




VaPro

Touch Free Hydrophilic Intermittent Catheter

Simple. Clean. Independence.



The VaPro and VaPro Plus hydrophilic intermittent catheters help provide the right balance of ease-of-use and protection, thanks to the unique design and exceptional features.



People First: Timm Seestadt

"I thought that racing wheelchair was amazing," Timm says. "It only had three wheels and looked really fast. I was very excited and I realized that's what I wanted to do – to race and to go as fast as I possibly could."



Patient Needs

Easy Handling

- Patients benefit from an intermittent catheter that can be easily gripped and manipulated, and the need is amplified for those with limited dexterity

Independence

- For patients, keeping or achieving their independence can be everything
- Balancing intermittent self-catheterization with daily life can be a challenge

Adherence to a Schedule

- Learning to catheterize can create apprehension and challenges^{1,2}
- Mastering the IC procedure and adhering to a catheterization schedule is important to helping patients maintain long-term health^{1,3}
- If the catheter is easy to use, the patient is more likely to adhere to a catheterization schedule³





What healthcare practitioners are saying about the VaPro catheter:⁴

“Easy to explain to patients how to use the product.”

— Medical Assistant

“Easy to maneuver and lubricity is a plus.”

— Registered Nurse


“Less steps means better adherence to proper use.”

— Registered Nurse



VaPro and VaPro Plus Catheters

- The VaPro catheter is designed to help reduce bacterial contamination and be easy to use
- In an independent assessment, 89% of respondents who evaluated the VaPro catheter were “likely” or “very likely” to recommend the VaPro catheter with “ease of use” and reduced contamination as the most often mentioned reasons⁴



For those who prefer gel lubrication

Advance Plus intermittent catheter system

The Advance Plus intermittent catheter system features a Touch Free technique enabling users to catheterize with confidence – virtually anytime, anywhere.



Available with an integrated 1000 mL collection bag

A4353 Defined⁵

An intermittent urinary catheter with insertion supplies (A4353) is a kit, which includes a catheter and all supplies necessary for a single, sterile insertion (see below).

Code A4353 may be used if any of the following 1, 2 or 3 is supplied:

1. A single sterile package containing both an intermittent urinary catheter and all necessary insertion/collection supplies; or,
2. A sterile intermittent urinary catheter plus a separately-packaged sterile kit containing all necessary insertion/collection supplies; or,
3. A sterile “no-touch” type of catheter system.

1 Achterberg, et al. Adherence to clean intermittent catheterization procedures: determinants explored. *Journal of Clinical Nursing* 2008;17:394-402

2 Logan K, Shaw C, Webber I, et al. Patients' experiences of learning clean intermittent self-catheterization: A qualitative study. *J Adv Nursing* 2008; 62:32-40

3 Shaw C, Logan K, Webber I, eal. Effect of clean intermittent self-catheterization on quality of life: A qualitative study. *J Adv Nursing* 2008; 61:641-50

4 Data on File. VaPro Sample Testing Healthcare Professionals Quantitative Research: Summary of Findings. Oct 2014

5 A4353 Source

Local Coverage Article for Urological Supplies - Policy Article - Effective December 2013 (A25230)

Article ID A25230

Article Title Urological Supplies - Policy Article - Effective December 2013

Contractor Name: NHIC, Corp.

Contract Type: DME MAC

DME MAC

Section: Coding guidelines

VaPro

Touch Free Hydrophilic Intermittent Catheter

VaPro Touch Free Hydrophilic Intermittent Catheter

Stock No	Length (inches)	Size	Box Qty
72082	8	8Fr	30
72102	8	10Fr	30
72122	8	12Fr	30
72142	8	14Fr	30
72084	16	8Fr	30
72104	16	10Fr	30
72124	16	12Fr	30
72144	16	14Fr	30
72164	16	16Fr	30
73124	16	12Fr Coudé	30
73144	16	14Fr Coudé	30
73164	16	16Fr Coudé	30

VaPro Plus Touch Free Hydrophilic Intermittent Catheter

Stock No	Length (inches)	Size	Box Qty
74122	8	12Fr	30
74142	8	14Fr	30
74124	16	12Fr	30
74144	16	14Fr	30

Advance Plus Kit Touch Free Intermittent Catheter System

Stock No	Length (inches)	Size	Box Qty
96064	16	6Fr	100
96084	16	8Fr	100
96104	16	10Fr	100
96124	16	12Fr	100
96144	16	14Fr	100
96164	16	16Fr	100
96184	16	18Fr	100
97124	16	12Fr Coudé	100
97144	16	14Fr Coudé	100
97164	16	16Fr Coudé	100

Advance Plus Pocket Touch Free Hydrophilic Intermittent Catheter

Stock No	Length (inches)	Size	Box Qty
93104	16	10Fr	100
93124	16	12Fr	100
93144	16	14Fr	100
93164	16	16Fr	100

- Ready-to-use, single-use intermittent catheter
- Hydrophilic-coated PVC catheter
- Protective tip and sleeve supports no-touch insertion technique
- 2 smooth catheter eyelets
- Color-coded funnel
- Adhesive tab for adhering packaging to a surface

This intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female, and pediatric patients who need to drain urine from the bladder.

WARNING: To help reduce the potential for infection and/or other complications, do not reuse. If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner.

PRECAUTION: Please consult a medical professional before using this product if any of the following conditions are present: Severed urethra, unexplained urethral bleeding, pronounced stricture, false passage, urethritis – inflammation of the urethra, prostatitis – inflammation of the prostate gland, epididymitis – inflammation of the epididymis (testicle tube). Self-catheterization should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterization frequency varies by person, the recommended frequency of your catheterization should be provided by your healthcare practitioner. For any other questions about your catheterization, please contact your healthcare practitioner.

CAUTION: Rx Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product. Prior to use of this device, be sure to read (i) the complete information on how to use this device including Warnings, Cautions, and Instructions for Use, and (ii) all other package inserts and labels supplied with the product.



Hollister Incorporated

2000 Hollister Drive
Libertyville, IL USA

For detailed clinical questions concerning our products:

1.888.740.8999
For orders only: 1.800.323.4060

www.hollister.com

The Hollister logo, Advance Plus, Advance Plus Pocket, VaPro, and VaPro Plus are trademarks of Hollister Incorporated.

©2014 Hollister Incorporated.
922678-1114